

Cameras and DVRs

Body Temperature Scanner Camera (Item: TMT2)

The TMT2 body temperature scanner will scan for above-normal temperatures in 1-2 seconds even when a person is wearing a mask. Each scan is saved in the cameras' software and the software can be integrated with a firm's network. By saving scans of each person, businesses can stay on top of worrisome symptoms among their employees and customers. Tracking scans can also contribute to contact tracing.



Features:

Get back to work safely. Give your customers or employees peace of mind with the touchless temperature check camera.

- Check all of your employees' body temperatures on arrival using facial recognition
- Keep records of the body temperature scans
- High precision infrared temperature detector
- Network to your existing access control system
- Can scan one person every 1-2 seconds
- Camera can stand-alone, no network needed
- Check all your customer's temperatures
- Support SDK and HTTP Protocols- **Coming soon**
- When mask is required, the camera automatically prompts you to put a mask on if attached to magnetic door lock will not allow entry without a mask.
- Build your employee data base or auto-register
- Email Alerts firmware v17.18.1.28
- 1 year warranty
- Optional stands available TMT2 Stand, TMT2 Stand2, TMT2 Stand3, TMT2TP

Perfect for nursing homes, schools, factories, warehouses, access gates, office buildings, hair salons, hospitals - any workplace or public area considering removing restrictions.

Returns and Refunds

To be eligible for a return, your item must be unused in its original unopened packaging and returned within 14 Days from date of sale. If

Cameras and DVRs

14 days have gone by since your purchase, unfortunately, we cannot offer you a refund or exchange. Due to health and safety precautions, used items CANNOT be returned for exchange or refund.

Exchanges (if applicable)

We only replace items if they are defective or damaged during shipping. If you need to return to exchange for the same item contact Customer Support. All returns are subject to final inspection upon receipt

This device is not intended for use in the diagnosis of any disease or other conditions or in the cure, mitigation, treatment, or prevention of any disease*

Body scanner video

FIRMWARE UPDATE

[v17.18.1.28](#) Click this link to open a new webpage. On the new webpage will be a **Download** button which when clicked will give you a choice to direct download the executable ZIP file named **17.18.1.28_VSTMT007.EXE** or to Save it to your Dropbox account.

Once you have downloaded the EXE file, double-click the file execute it, you will be prompted for a location to unzip the contents of the file, default location is "c:TMT2_Firmware_Update", click the Unzip button to continue. After you have extracted the file contents close the zip window then browse to folder location and click to open the folder then open the folder 17.18.1.28_VSTMT007, in that folder there is a file named Upgrade.Directions.txt that explains how to install the firmware update.

Update includes various bug fixes, all new English voice prompts

This device is not intended for use in the diagnosis of any disease or other conditions or in the cure, mitigation, treatment, or prevention of any disease*

Cameras and DVRs

Tech Specs:

- Power consumption 20W(MAX)
- Power input DC 12V/3A
- Interface connections:
 - USB port 1USB port (currently not supported) coming soon
 - Alarm output 1 channel relay output
 - Weigand port Support input/output 26 and 34
 - Internet interface RJ45 10M/100M Ethernet
- Temperature:
 - Range 86F- 113F
 - Response time:.5 seconds
 - Best distance: 1 1/2 feet
 - Accuracy .3 +or -
- Camera:
 - Color Camera .01 LUX
 - 2Mega pixel 1920*1080
 - Color Camera .01 LUX
 - WDR
 - Camera: Support live detection
 - Face Recognition:
 - Height 1.2-2.2 M, angle adjustable
 - View angle Vertical ±30 degree
 - Support 22400 faces database and 100000 records
 - Distance 2-10 Ft.

Includes:

- Pole mount
- Power adaptor